

Central sterile services department (CSSD)



H2O2 Low Temperature Plasma Sterilizer (Dual Chamber)

200L



Closed-end Traceability

Connect to rapid readout biological reader for closed-end traceability of sterilization results and biological incubation results

Smart Recognition

Equipped with RFID reader, camera and bar code reader to realize fast recognition of consumables and auto cycle binding. Assisting users to fast login according to access level settings.

Efficiency Improvement

Shortest cycle time 20mins.

H2O2 low temperature plasma sterilizer

100L



AI Processor

Compatible with MS Windows® and Android® system for remote OTA system update. Voice interactivity for voice control of routine operations.

Multiple Cycles

6 preset cycles can be run simultaneously in the dual chambers

Dual Chamber

Independent dual chamber design to run two different cycles simultaneously to drastically improve turnover of surgical devices.



Automatic washer disinfector

Fully automatic surgical instrument cleaning and disinfection equipment, fully automated processing of cleaning, disinfection, drying, and oiling:1. Volume: 500L2. Material: 304 stainless steel3. Double door interlocking, automatic lifting door4. Water consumption: 23L/step5. Cleaning time: 27 minutes6. Meet the B-level disinfection standard of German BGA7. Working mode: electric heating/steam heating/dual energy of steam and electricity8. Configuration: 1 host, 1

pulse vacuum steam sterilizer

CSSD must be equipped with high-temperature sterilization equipment:1. Capacity: 600L2. Door: Double door pressure safety interlock device, motorized door opening method3. The main inner shell, reinforcement ribs, and door panels are all made of 304 stainless steel; Design life of 15 years4. Working mode: Electric heating5. Configuration: 1 host, 1 disinfection vehicle, and 2 transport vehicles6. Optional: 1 air compressor



H2O2 Low Temperature Plasma Sterilizer

4 preset cycles: Short, Standard, Endo, Drying
Advanced features: H2O2 detection, soft heating, upgraded vacuum
Premium build: Aluminum chamber, retractable shelves, 10" touch screen
Quiet & safe: ≤55dB, protects sensitive instruments
Connectivity: Compatible with HIS & CSSD systems



info@dijla-medical.com
www.dijla-medical.com



Erbil / italy2 / office279
Baghdad / ALYarmouk



07755459000
07755459000