

**A TRUE HYBRID.
CART-BASED PERFORMANCE;
HAND-CARRIED DEVICE.**



TORSO-ONE™

A TRUE HYBRID.

**CART-BASED PERFORMANCE;
HAND-CARRIED DEVICE.**

EchoNous™ developed the Kosmos Torso-One to give point-of-care clinicians the most information possible at the patient's bedside. By infusing the power and performance of premium ultrasound systems, combined with industry-leading AI, into a sophisticated hand-carried form factor, we have achieved that goal - though we've only just begun.

Built for real-life POCUS situations

Get optimal imaging without lifting a finger. Kosmos, with its superior image quality and Auto Preset which adjusts the exam type in real-time based on the anatomy you are scanning — you can spend less time trying to optimize the image and more time caring for your patients.

CW, PW, TDI and Auto Doppler

With CW, PW, and TDI capabilities, you can acquire valuable, in-depth information at the bedside to make more confident patient management decisions. Auto Doppler goes even further! Available with PW and TDI, it will place the sample gate in the optimal location for more efficient and consistent cardiac valve interrogation.

One Probe, Multiple Configurations



TORSO-ONE SPECS

- **AI Features:** Trio, Auto EE, AI FAST, Auto
- **Doppler, Auto Preset, AI VTI Trace**
- **Exam Types:** Cardiac, Lung, & Abdomen
- **Modes:** 2D, M-Mode, Color, PW, CW, TDI
- **Doppler Measurements:** Comprehensive
- **Doppler calculations**

- **2D Frequency:** 1.5 - 4.5 MHz
- **Center Frequency:** 3 MHz
- **Footprint:** 20.7mm x 16mm
- **Mass:** 275g
- **Cable Length:** 1.5 m



info@dijla-medical.com
www.dijla-medical.com



Erbil / italy2 / office279
Baghdad / ALYarmouk



07755459000
07755459000